LOOP BAND INSTRUCTION MANUAL

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Safety Instructions – Please Note

- Natural product with reliable linear resistance behaviour. Sanctband | Sanctband Active is a functional piece of exercise device for physical therapy, muscle increase and fitness.
- ATTENTION: Sanctband | Sanctband Active are made from natural latex which can cause allergic reactions.**
- When using Sanctband | Sanctband Active, always follow the recommendations of your doctor or therapist.
- Children may use Sanctband | Sanctband Active only under adult supervision and after a therapy-accompanying introduction. Not suitable for infants.
- Check Sanctband | Sanctband Active for damages (holes, cracks, cuts) prior to every use. Do not use damaged loop band, replace them instead.
- Do not overstretch Sanctband | Sanctband Active ($F \geq 200\%$), do not create too much tension. Choose the loop band that matches your physical condition and muscle strength.
- Be careful with the exercises at head height! Keep away from eyes, the loop band could slip out of hand or holder and snap in the face!
- Ensure the loop band cannot release from door frame or loosen from the chair or other objects. Ensure that door is shut while exercising.
- Sanctband's | Sanctband Active's Recommended Cleaning Method: Clean with wet clean cloth.
- Disposal of the devices: Follow local disposal procedure as per user facility.**

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