TUBING WITH HANDLES

INSTRUCTION MANUAL

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LATIHAN KETAHANAN

Safety Instructions – Please Note

1. Natural product with reliable linear resistance behaviour. Sanctband | Sanctband Active is a functional piece of exercise device for physical therapy, muscle increase and fitness.
- 2. ATTENTION: Sanctband | Sanctband Active are made from natural latex which can cause allergic reactions.**
3. When using Sanctband | Sanctband Active, always follow the recommendations of your doctor or therapist.
4. Children may use Sanctband | Sanctband Active only under adult supervision and after a therapy-accompanying introduction. Not suitable for infants.
5. Check Sanctband | Sanctband Active for damages (holes, cracks, cuts) prior to every use. Do not use damaged tubings, replace them instead.
6. Do not overstretch Sanctband | Sanctband Active ($F \geq 200\%$), do not create too much tension. Choose the tubing that matches your physical condition and muscle strength.
7. Be careful with the exercises at head height! Keep away from eyes, the tubing could slip out of hand or holder and snap in the face!
8. Ensure the tubing cannot release from door frame or loosen from the chair or other objects. Ensure that door is shut while exercising.
9. Sanctband's | Sanctband Active's Recommended Cleaning Method: Clean with wet clean cloth.
- 10. Disposal of the devices: Follow local disposal procedure as per user**

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