

Figure: Force-Elongation of Sanctband Exercise Tubing - ID 5.0mm




Table: Force-Elongation for Sanctband Exercise Tubing - ID 5.0mm (Force in Pounds)

Elongation %	Peach	Orange	Lime Green	Blueberry	Plum	Gray
25%	1.61	1.87	2.71	3.19	3.41	5.02
50%	2.66	3.12	4.55	5.32	5.94	8.29
75%	3.43	4.05	5.83	6.66	7.74	10.52
100%	4.05	4.75	6.91	7.92	9.11	12.52
125%	4.58	5.32	7.79	8.98	10.32	14.26
150%	5.06	5.90	8.60	9.99	11.57	15.71
175%	5.57	6.38	9.31	10.89	12.63	17.01
200%	6.03	6.91	10.08	11.79	13.73	18.28
225%	6.53	7.44	10.71	12.78	14.89	19.76
250%	6.95	7.90	11.44	13.49	16.10	21.16
275%	7.39	8.29	12.06	14.30	17.03	22.31
300%	7.92	8.87	12.87	15.25	18.28	23.72