

Figure: Force-Elongation of Sanctband Exercise Tubing - ID 6.5mm




Table: Force-Elongation for Sanctband Exercise Tubing - ID 6.5mm (Force in Pounds)

Elongation %	Peach	Orange	Lime Green	Blueberry	Plum	Gray
25%	2.11	3.10	3.63	4.20	4.53	5.39
50%	3.65	5.46	6.14	7.02	7.59	8.84
75%	4.69	6.78	7.83	8.76	9.79	11.29
100%	5.61	8.14	9.35	10.34	11.75	13.62
125%	6.34	9.35	10.56	11.75	13.29	15.49
150%	7.04	10.34	11.66	12.96	14.89	17.05
175%	7.72	11.44	12.78	14.08	16.48	18.85
200%	8.36	12.39	13.79	15.16	17.91	20.39
225%	9.09	13.51	15.00	16.61	19.69	22.44
250%	9.68	14.54	15.97	17.49	21.03	23.80
275%	10.27	15.55	16.98	18.59	22.53	25.43
300%	11.07	16.63	18.26	19.73	24.29	27.24