

Figure: Force-Elongation of Sanctband Loop Bands




Table: Force-Elongation for Sanctband Loop Bands (Force in Pounds)

Elongation %	Peach	Orange	Lime Green	Blueberry
25%	1.34	2.27	3.04	4.22
50%	2.35	4.03	5.20	7.24
75%	3.01	5.21	7.04	9.70
100%	3.70	6.29	8.36	11.55
125%	4.25	7.26	9.70	13.42
150%	4.82	8.14	10.91	15.20
175%	5.28	9.00	12.12	16.94
200%	5.70	9.86	13.27	18.68
225%	6.18	10.76	14.52	20.66
250%	6.66	11.62	15.75	21.89
275%	7.06	12.39	16.87	23.43
300%	7.50	13.27	18.24	24.60